

# Sample Bottle Assembly



## One of a Kind

Aquasyn's unique sample bottle assembly with an integrated Valbow® solves a number of process problems and has become an industry standard. It is designed for those that prefer a reusable component over one that is disposable.

Sample bottles can be provided with any number of options including different bottle sizes, color coded caps, and unusual Valbow® sample valve configurations.

Unlike other designs, Aquasyn also offers inspection and refurbishment services for our sample bottle assemblies.

## Features

- Bottle sizes from 50ml to 1000ml
- Color coded caps
- Fully autoclavable
- Custom designs available
- Inventory of stocked product
- Industry standard acceptance
- Lightweight design

## Valbow®



## Unique Products

The Valbow® was created some years ago as a solution to one of our customer's unique problems. It has evolved into a popular Aquasyn product and is stocked in a variety of sizes. The Valbow® is uniquely suited for that one spot where nothing else will do. Whether it's a hard to reach sample valve on an R/O system, or for a directional change valve on a Bioreactor, valbows just seem to come in handy.

## Features

- Unique problem solving design, unlike any other in the industry
- 90 degree angle pattern for directional changes
- Minimal hold up volume
- Sampling capability
- Hygienic clamp, orbital weld, end cap and other connections available
- Custom designs and other special options available

