



Ordering Information

Valve to Valve Combinations




SAL1 08 L 2 08 09 09 09 E 19 E 16 Sterile Access Valve to Valve Combination



- Type: Sterile Access Valve Left Side (SAL1)
- 1st Valve: 1/2"
- Material: 316L Stainless Steel
- Finish: 15Ra/Electropolished
- 2nd Valve: 1/2"
- Outlet A: 1/2" Tubing
- Outlet B: 1/2" Tubing
- Outlet C: 1/2" Tubing
- Diaphragm 1: EPDM
- Operation 1: Actuator Normally Closed (N/C)
- Diaphragm 2: EPDM
- Operation 2: SS Manual Bonnet Assembly

HGV1 16 L 7 12 AE AE AD P 16 P 16 GMP Valve to Valve Combination



- Type: Horizontal GMP Valve (HGV1)
- 1st Valve: 1"
- Material: 316L Stainless Steel
- Finish: 20Ra (No Electropolish)
- 2nd Valve: 3/4"
- Outlet A: 1" Sanitary Clamp
- Outlet B: 1" Sanitary Clamp
- Outlet C: 3/4" Sanitary Clamp
- Diaphragm 1: PTFE
- Operation 1: SS Manual Bonnet Assembly
- Diaphragm 2: PTFE
- Operation 2: SS Manual Bonnet Assembly

¹Hastelloy is a trademark of Haynes International
²AL6XN is a trademark of Allegheeny Ludlum

Selection Chart

