



Key of Table

No visible effect	
No data	
Cracked or crazed	

Table 5 **Desinfectants**

Trade Name	Manufacturer	Temp. in Celsius	Concentration in Celsius	Time in Hr	Applied Strain								
					0.0	0.2	0.4	0.6	0.8	1.0	1.2	1.4	
Aldetex R -28	ICI Pharma	23	100	72									
AseptisolR	Dr. Bode & Co.	23	2.5	72									
Bodephen R	Dr. Bode & Co.	23	3	72									
BomixR	Dr. Bode & Co.	23	3	72									
CIDEXR	Johnson & Johnson	23	100	72									
DettolR	Reckitt & Colemann	23	5	72									
Gigasept R FF	Schulke & Mayer	23	15	72									
GrotanatR	Schulke & Mayer	23	4	72									
H.A.C.	ICI Pharma	23	1	72									
Habitane R	ICI Pharma	23	1	72									
Harvey's Vapo-Steril R	Barnstead/Thermol	135	100	6.6									
Lysetol R FF	Schulke & Mayer	23	4	72									
Neodish er R A8	Dr Weigert	80	5	72									
S&M Labor	Schulke & Mayer	23	1	72									
Sagrotan R MED	Schulke & Mayer	23	0.5	72									
Secumatic R FHZ	Henkel	60	0.3	144									
Secumatic R FRE	Henkel	60	0.5	144									
Secumatic R FORTE	Henkel	23	0.5	72									
SokrenaR B		60	2.5	72									